

# Clinical Record Requirements for 3M™ Filtek™ Matrix

This guide will help you verify the criteria for getting your case approved and avoid delays in delivery of the custom matrix.

## Scan Requirements



Full arch scans that reflect current state of dentition.

- ① Upper Arch
- ② Lower Arch
- ③ Right Bite
- ④ Left Bite

## Criteria for Acceptance or Rejection

### Required

- ☑ All necessary information captured on scans
- ☑ All incisal edges and marginal ridges
- ☑ Clear gingival margin with 3 to 5 mm of lingual and buccal gingiva

### Optional

- ⊖ Aligner attachments are acceptable within scans but will be digitally removed and must be removed before seating matrix
- ⊖ Outside lab wax up scans can be included to use as reference for the restorative design

**Note:** To ensure reliable fit and function of the Matrix, outside lab wax ups cannot be precisely replicated and may have to be adjusted based on Filtek Matrix requirements.

### Rejected

- ⊘ Scans older than 90 days
- ⊘ Outside lab wax up scans only without other records
- ⊘ Scans containing partial dentures on restorative arch
- ⊘ Double surfacing on scan
- ⊘ Incomplete capture of margins
- ⊘ Incomplete capture of the incisal edge(s)
- ⊘ Brackets or bonded lingual retainers present on restorative arch
- ⊘ Missing and/or excess data

## Required Photos

All photos are captured and are clear



① Right Lateral

③ Upper Occlusal


⑤ Frontal

② Left Lateral

④ Lower Occlusal


## Optional Photos

Full facial smile photo can be included to help set orientation to ensure proper leveling



⑥ Lip at Rest (optional)

⑦ Full Facial Smile (optional)



**We're here for you.**

Call the 3M Oral Care Customer Care team for troubleshooting and technical support: **1-800-634-2249**.

**3M**

3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA  
Phone 1-800-634-2249  
Web 3M.com/dental

3M and Filtek are trademarks of 3M. © 2023, 3M.  
All rights reserved. 70-2013-8155-8