

# 3M<sup>™</sup> Attest<sup>™</sup> Auto-reader 490 System<sup>\*</sup> Quick Reference Guide

# Quick Reference Guide



#### Inspect

Remove BI from pouch and reseal. Verify media ampoule is intact and indicator stripes on cap are blue.



#### **Pouch**

Place BI and chemical indicator (CI) in peel pouch indicated for use in vaporized hydrogen peroxide (VH2O2) sterilization processes.



#### **Process**

Seal pouch and place in sterilizer with white side facing up. Process load as recommended by sterilizer manufacturer.



### Identify

After completion of cycle, inspect and verify media ampoule is intact and indicator stripes on cap have turned towards pink. Identify the BI by writing on the label.

# STEAM and VH2O2 BI Activation and Incubation



Crush

Squeeze BI in activator to close cap and crush media ampoule.



Flick

Remove BI and flick wrist to ensure media flows to growth chamber.



Inspect

Verify media has transferred to BI growth chamber at the vial bottom.



#### Incubate

Place activated BI in any well. Display shows remaining minutes of incubation or test results. DO NOT REMOVE BI until the (+) or (-) symbol indicates the test is complete.

# 3M<sup>™</sup> Attest<sup>™</sup> Auto-reader 490

#### Positive Control BI


For both Steam and VH2O2, each day a processed BI is incubated, activate and incubate an unprocessed BI as a positive control. The control BI should be from the same lot code as the processed BI.

#### Results After Incubation

Control BI: Positive (+) Acceptable

Processed BI: Negative (-) Acceptable Sterilization
Processed BI: Positive (+) Failed Sterilization


Remaining Incubation Minutes





#### NOTES:

- Activate and incubate 1295 BI within 1 hour of completion of the VH2O2 sterilization cycle.
- Gloves are recommended when activating 1295 Bl.
- Allow 1491 and 1492V Bls to cool 10 minutes prior to activation.

## **Interpretation of Solventum Sterilization Assurance Core Products**







STEAM Indicator

3M™ Attest™ eBowie-Dick
Test System

3M™ Attest™ eBowie-Dick
Test Card (left)

3M™ Attest™ eBowie-Dick
Auto-reader (right)

Unprocessed

36. Attest

data Cale to Cale

and Cale to Cale

and Cale to Cale

and Cale

Processed

SM Attest

Salar Attent

Salar At

circle on test card darkens





Solventum Corporation 3M Center, Building 275 2510 Conway Avenue East Maplewood, MN 55144 © Solventum 2024. Solventum, the S logo are trademarks of Solventum or its affiliates. 3M, the 3M logo, Attest and Comply are trademarks of 3M. Other trademarks are the property of their respective owners. 70-2011-7096-9 Rev 5/24

\*Refer to 3M<sup>™</sup> Attest<sup>™</sup> Auto-reader 490 Operator Manuals and 3M<sup>™</sup> Attest<sup>™</sup> Biological Indicator 1295, 1491 and 1492V Instructions for Use for complete instructions and other important information. Refer to the indications for use should align with the sterilization cycle to be monitored. This Quick Reference Guide applies to 3M<sup>™</sup> Attest<sup>™</sup> Auto-reader 490 having software version 4.2.7 or greater. Refer to the Instructions for Use for complete instructions and other important information related to the use of all other products included in this Guide.



