

Enabling Process Simplification through Chromatographic Clarification with High Titer Cell Culture

David Chau¹, Aaron Dryden², Mark Smith², and Matt Zustiak³.


¹3M Center, Maplewood, MN, USA. 55114. ²Thermo Fisher Scientific. 1325 North 1000 West. Logan, UT, USA. 84321.

³Thermo Fisher Scientific. 4766 LaGuardia Drive. St. Louis, MO, USA. 63134.

Introduction


With advancements in upstream bioprocesses resulting in higher titer and cell densities, newer technologies has had to meet the demand for mass-transfer and mixing needs of an intensified process. However, as a new framework has been developed to support intensified cell cultures, newer challenges have emerged for the purification stages. The combination of higher cell levels along with process-related impurities (e.g., DNA/HCP/etc.) have posed a considerable challenge for downstream operations to meet the advancements that have been made in upstream bioprocesses.

Figure 1. Cell Culture Contaminants and Process-Related Impurities Plotted Along Their Size and Isoelectric Point.



3M™ Harvest RC is a harvest separation solution that utilizes fibrous anion exchange (AEX) chromatography to efficiently separate the cells, cell debris, and DNA from the harvest fluid containing the target product¹. 3M™ Harvest RC is a new single-stage, single-use chromatographic clarification solution, designed as an efficient option to increase monoclonal antibody (mAb) yields while streamlining the upstream process by replacing the centrifuge and/or depth filtration process steps.

Figure 2. Process Compression using 3M™ Harvest RC Chromatographic Clarifier




Results

Performance of 3M™ Harvest RC Chromatographic Clarifier vs. Traditional Depth Filtration

A cell culture was grown and harvested at different days to illustrate the impact of viability and pack cell volume on the 3M™ Harvest RC Chromatographic Clarifier and a single stage depth filter technology. The cell culture was harvested on day 7, 9 and 12 to illustrate the impact of viability on the performance of clarification. As viability dropped over the course of the culture run, material was harvested and challenged onto the two different clarification technology.

Figure 3. Throughput versus Pack Cell Volume for 3M™ Harvest RC Chromatographic Clarifier Performance compared to Traditional Depth Filtration Technology.



Robust Scalability of 3M™ Harvest RC Chromatographic Clarifier

Figure 4. Throughput Performance of 3M™ Harvest RC Across Different Scales





Figure 5. Turbidity Reduction using 3M™ Harvest RC Across Different Scales



Case Studies with Different Cell Cultures to Robust Performance

Different cell cultures were grown by Thermo Fisher Scientific to challenge the 3M™ Harvest RC Chromatographic Clarifier technology under a broad range of culture conditions. These cultures had diverse characteristics in cell density, viability, and pack cell volume. These different cultures were used to challenge the 3M™ Harvest RC as well as a traditional dual stage depth filter train from Supplier A. The results were then used to calculate the total surface area that would be needed to conduct a 2000L and 5000L scale run and the materials that would be needed.

Table 1. Four Different Cell Culture Characteristics for Challenging 3M™ Harvest RC Chromatographic Clarifier.

Culture	1	2	3	4
VCD (x 10 ⁶ cells/mL)	27	20	15	9
Viability (%)	91	86	83	76
PCV (%)	12	14	8	4

Table 2. Total Surface Area Needed for Clarification at 2000 and 5000L Run Using A Traditional Dual Stage Depth Filtration Process versus 3M™ Harvest RC Chromatographic Clarifier.

Culture	1	2	3	4
Dual Stage Depth Filter	50 (total)	57 (total)	27 (total)	30 (total)
3M Harvest RC	29	41	19	14
5K filter area (m ²)	123 (total)	142 (total)	66 (total)	74 (total)
3M Harvest RC	74	102	49	35

Table 3. Case Study for Culture #3 for a 5000L Scale Run.

	Supplier A		3M™ Harvest RC	Improvement with 3M™ Harvest RC	
	Stage 1	Stage 2	Single stage	Single-stage reduction	
EFA (m ²)	44	22	49*	26% reduction	
No. of capsules	40	20	45	25% reduction	
Footprint	-	-	11.4	7.6	33% reduction
Pre-use flush (L)	4,400 (100 L/m ²)	2,200 (100 L/m ²)	6,600 L	1,225 L (25 L/m ²)	81% reduction
Yield (%)	-	86	100	14% increase	

Conclusion

- The 3M™ Harvest RC Chromatographic Clarifier can act as a single stage single use technology for promoting process compression through an innovative chromatographic clarification approach for capturing cells, cell debris, and soluble impurities.
- The results illustrated that the 3M™ Harvest RC had superior capacity compared to a traditional single stage or dual stage depth filter technology which can lead to reduced set-up time and waste generation.
- The use of the 3M™ Harvest RC could result in highly clean filtrate with a single stage that would result in a simplified process development and process set-up. This single stage set-up also resulted in improved yield compared to traditional multiple stage set-up that would translate to product loss.

References

- Almeida, A., Chau, D., Coolidge, T., El-Sabbahy, H., Hager, S., Jose, K., Nakamura, M. and Voloshin, A., 2021. Chromatographic capture of cells to achieve single stage clarification in recombinant protein purification. *Biotechnology Progress*.

Acknowledgement

We would like to thank Alexei Voloshin and Masa Nakamura for their support and guidance during the planning and experimentation stages. We would also like to thank Qingwei Luo, Patrick Bennet, Michelle Nolasco, Kristina Pleitt, and Jon Kruger from Thermo Fisher Scientific for their support and work in providing the material and experimental work.