

## 3M™ Nasogastric Securement Device


Designed to deliver reliable securement to help keep nasogastric tubes in place with an easy-to-use solution that helps support pressure injury prevention practices and minimise the risk for medical adhesive-related skin injury (MARS).<sup>1,2</sup>

### Non-adhesive section at the base of the nose

Intended to help the tube hang freely to reduce contact between it and patient skin.

### Transparent for visibility

Assess skin for complications and monitor tube positioning without removal.



### Fibre mesh backing

Provides extra strength.

### Non-adhesive tab for removal

Easy removal  
in one piece.



**Skin performance:** The feather designates products that deliver the securement power you need while minimising damage to skin.

### Reliable securement

Unique adhesive secures nasogastric tubes to dry skin and lasts for up to 48 hours with minimal migration even after becoming wet and oily.

### Ease of use

Two-piece liner allows for easy application with gloves; non-adhesive tab at the end of the tail allows for easy removal. Supports standardisation of technique.

### Gentle to skin

Non-sensitising adhesive has low potential to cause dermatitis and minimises pain and damage to skin upon removal.<sup>1,2</sup>

### Designed to help reduce pressure injury risk

Non-adhesive area is designed to minimise pressure at the base of the nares, while product transparency facilitates skin assessment.


1. 3M data on file. EM-05-014947.  
2. 3M data on file. EM-05-014954.

The securement solution you choose for your patients – whether a device or tape – should keep tubes in place while minimising migration to help avoid disrupting nasogastric (NG) therapy. Consider how the 3M™ Nasogastric Securement Device performed in a comparative study.

## A comparison of the 3M™ Nasogastric Securement Device to the Hollister® Feeding Tube Attachment Device.<sup>1</sup>

### More reliable securement

3M™ Nasogastric Securement Device had a statistically higher percentage of devices remaining on the skin than Hollister® Feeding Tube Attachment Device after 48 hours and 72 hours of wear.\*



### Less migration

The 3M™ Nasogastric Securement Device showed 60% less migration than the Hollister® Feeding Tube Attachment Device at 72 hours.\*†





The secure choice. Trust the 3M™ Nasogastric Securement Device to offer consistent performance that is strong, yet gentle.

To order, contact your 3M Health Care Sales Representative or call the 3M Health Care Customer Helpline at 0870 536 0036. Visit [3M.co.uk/Medical](http://3M.co.uk/Medical) for more information.

Refer to product Instructions for Use (IFU) for other important information.

### Ordering information

Product	Product name	Product number	Size	Devices/box	Boxes/base
	3M™ Nasogastric Securement Device, Large	1501NG	5.7cm x 7.3cm	50	4
	3M™ Nasogastric Securement Device, Small	1500NG	4.5cm x 4.7cm	50	4

#### Medical Care Solutions Division

3M United Kingdom PLC Charwood Campus 10 Bakewell Road Loughborough LE11 5RB +44 (0)1509 611 611	3M Ireland The Iveagh Building Carrickmines Park Carrickmines Dublin 18 00 353 (01) 280 3555
---	---

\*As tested in a comparative, prospective clinical study of 64 healthy human subjects over 3 days of wear.

†A 50-gram weight was hung from each device for 15 minutes daily to simulate tube presence. Percentages were determined by the difference between devices with migration 30% (3M) and 78% (Hollister).