

3M™ Ranger™ Pressure Infusor System
High Flow Solutions



Next Generation Pressure Infusion

Broader fluid capacity. Advanced alarm system.

For use with 3M™ Ranger™ blood and fluid warming system.
Provides pressure to enable high flow fluid volumes up to
30,000 mL/hr.

Lead
The Way

3M

3M™ Ranger™ Pressure Infusor System

The 3M™ Ranger™ pressure infusor system is designed for use with the 3M™ Ranger™ blood/fluid warming unit and 3M™ Ranger™ high flow warming sets. The pressure infusor system can accept 250 mL to 1000 mL fluid bags and provides maximum dynamic operating pressure of 300 mmHg.

The pressure infusor system includes three separate pieces: the pressure infusor, IV pole and IV pole base.



3M™ Ranger™ Pressure Infusor

Model 14500

For Use With:	Pressure infusor IV pole (model 90068) and Ranger IV pole base (model 90124)
Dimensions:	40 L x 51 W x 20 H cm (15.75 L x 20 W x 7.75 H in)
Weight:	7.7 kg (17 lbs)
Case Quantity:	1/Box

3M™ Ranger™ Pressure Infusor IV Pole

Model 90068

For Use With:	Ranger pressure infusor (model 14500) and Ranger pressure infusor base (model 90124)
Dimensions:	187 L x 11 W x 4 H cm (73-1/2 L x 4-1/2 W x 1-3/4 H in)
Weight:	2.8 kg (6.2 lb)
Case Quantity:	1/Box

3M™ Ranger™ Pressure Infusor IV Pole Base

Model 90124

For Use With:	Ranger pressure infusor (model 14500) and Ranger pressure infusor IV pole (model 90068)
Dimensions:	81 L x 76 W x 20 H cm (32 L x 30 W x 8 H in)
Weight:	16.1 kg (35.5 lb)
Case Quantity:	1/Box



Infection Prevention Division 3M Health Care

3M Center, Building 275-4E-01
St. Paul, MN 55144-1000
U.S.A.
800-228-3957
3m.com/infectionprevention

3M is a trademark of 3M Company, used under license in Canada.
BAIR HUGGER, RANGER, and the BAIR HUGGER and RANGER logos are trademarks of 3M, used under license in Canada. Please recycle. © 3M 2014. All rights reserved. 602305D 2/14