
























# ECG Electrodes Product Comparison Chart

## Red Dot™ Electrodes

### The Right Electrode Matters

Choosing the right electrode can be difficult. That's why we've optimized our 3M™ Red Dot™ ECG Electrode portfolio to feature the options that meet your most common patient and performance needs. With 3M, optimal adhesion for a variety of patients, consistent and reliable trace quality, and ease of use come standard—helping you make a positive impact on patient care. **Reach for the Red. ●**



	Adult				Adult/Pediatric		Neonatal		
	Universal					Diagnostic			
	 <b>2228</b>	 <b>2560</b>	 <b>2660</b>	 <b>2239</b>	 <b>2268</b>	 <b>2330</b>	 <b>2269T</b>	 <b>2282</b>	 <b>2283</b>
	● Diaphoretic	● Diaphoretic	● Repositionable ● Fragile Skin (2660 only)*	● Fragile Skin*	● Fragile Skin*	● Fragile Skin (2330 only)*			
<b>Adhesion</b>			2660:  2670: 		 ● Gentle	2330:  2360: 			
<b>Trace Quality</b>	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●	●●●●●●●●
<b>MR Conditional</b>	2228: No 2244: Yes	2560: No 2570: Yes	Yes	No	Yes	Yes	Yes	No	No
<b>Radiolucent</b>	2228: No 2244: Yes	2560: No 2570: Yes	Yes	No	Yes	Yes	Yes	Yes	Yes
<b>Backing</b>	Foam	Foam	Soft Cloth	3M™ Micropore™ Tape	Soft Cloth	Polyester Film	Clear Tape	Soft Cloth	Soft Cloth
<b>Bag Quantity</b>	50	3, 5, 50	3, 5	50	3, 5	100	3	3	3
<b>Opened Pouch Expiration</b>	30 days	45 days	Use immediately	30 days	Use immediately	7 days	Use immediately	Use immediately	Use immediately