

3M™ Comply™ Lead Free Steam Indicator Tape

Visual assurance for maximum sterilisation.

Exposure monitoring products allows anyone to see at a glance whether their packs have been exposed to the sterilisation process.



Introducing a reliable 55m steam indicator tape that provides immediate information if the pack has been exposed to adequate steam sterilisation.

3M™ Comply™ Indicator Tape 1322

Disposable and reusable wraps.

3M™ Comply™ Indicator Tape 1355


Disposable wraps only.



Both provide a high level of adhesion – 1355 proving significantly high post-sterilisation adhesion to wrap.

Reliable steam indicator tape's that provides immediate information if the pack has been exposed to adequate steam sterilisation.

Process once exposed to steam.



The chemical indicator lines will turn dark when exposed to a steam sterilisation process.

Visual assurance for maximum sterilisation



Lead free



Over 50m in length
55m.



Convenient

Easy to use, single packaging.



No added latex



Reliable adhesives

Adheres to packaging material.



Easily identifiable

Dark black colour change and changes at the correct point in the cycle.

Description	Colour	Code	Packaging	Length
Lead Free Steam Indicator Tape Disposable and reusable wrap	Tan	1322 – 18mm	28 rolls/case	18mm x 55m
		1322 – 24mm	20 rolls/case	24mm x 55m
Lead Free Steam Indicator Tape Disposable wraps	Blue	1355 – 18mm	28 rolls/case	18mm x 55m
		1355 – 24mm	20 rolls/case	24mm x 55m



Scan or click on the QR code to watch the adhesion comparison demonstration.

**For more information call customer services
on 08705 360036 or visit [3M.co.uk/medical](https://www.3m.co.uk/medical)**