



CuROS™

Disinfecting Port
Protectors

Cap it with the right colour.

Use this CuROS Cap:

Only on:


Do not use on:



3M™ CuROS™ Disinfecting Cap
for needleless connectors



Needleless
connectors



Open female luers
(stopcocks or catheter hubs)



3M™ CuROS™ Stopper Disinfecting Cap
for open female luers



Open female luers
(stopcocks or catheter hubs)




Needleless
connectors



3M™ CuROS™ Disinfecting Cap
for Tego® Hemodialysis Connectors



Tego® Hemodialysis Connectors



Open female luers
(stopcocks or catheter hubs)



3M™ CuROS Tips™ Disinfecting Cap
for male luers



Male luers