



3M™ Ioban™ 2 Antimicrobial Incise Drape Selection Guide

Choose a powerful barrier against contamination.

Selection guide

Choosing a drape for your application.

Use the charts on these pages to select the best drape for common procedures. In some cases there may be more than one 3M™ loban™ 2 Antimicrobial Incise Drape that would work for your application. If this is the case, see the chart on the back page to compare specialty and general drape options for each surgery.













Easy-delivery

(products denoted with EZE) features include:

1. Polyethylene liner comes off smoothly, effortlessly, with no tearing, for easy drape delivery.
2. Drape and liner feature a full-width handle for control of liner release and drape application.
3. Strip at one end imprinted with stop signs indicates when the end of the drape is near.

	3M code	Description	Image	Size*	This size product may be useful for procedures such as:
Orthopaedics	6635	Incise drape		10cm x 20cm	Small limb, fingers, toes, paediatrics
	6661EZE	Framed delivery incise drape		26cm x 20cm	Ankle, wrist, paediatrics (can be applied by single user)
	6640EU	Small incise drape		35cm x 35cm	Small limb, ankle, wrist, spinal surgeries, paediatrics
	6640EZE	Small incise drape		35cm x 35cm	Small limb, ankle, wrist, spinal surgeries, paediatrics
	6650EU	Medium incise drape		56cm x 45cm	Knees, shoulders, hips, paediatrics
	6650EZE	Medium incise drape		60cm x 45cm	Knees, shoulders, hips, paediatrics
	6648EU	Large incise drape		56cm x 60cm	Knees, shoulders
	6648EZE	Large incise drape		60cm x 60cm	Knees, shoulders
	6651EU	Extra-large incise drape		56cm x 85cm	Hips, entire extremity/leg
	6651EZE	Extra-large incise drape		60cm x 85cm	Hips, entire extremity/leg
	6617	Isolation drape with loban 2 incise film and pouch		320cm x 213cm	Fractured neck of femur, hemi-arthoplasty, total hip replacement, intramedullary nailing
	6619	Large isolation drape with loban 2 incise film and pouch		328cm x 254cm	Fractured neck of femur, hemi-arthoplasty, total hip replacement, intramedullary nailing

*All dimensions are approximate.

	3M code	Description	Image	Size*	This size product may be useful for procedures such as:
Obstetrics and gynecology	6648	Large incise drape		56cm x 60cm	Open abdominal/pelvic surgery
	6648EZE	Large incise drape		60cm x 60cm	Open abdominal/pelvic surgery
	6657	C-section pouch with loban 2 incise film		87cm x 74cm	Caesarean section
Neurology	6640EU	Small incise drape		35cm x 35cm	VP shunt
	6640EZE	Small incise drape		35cm x 35cm	VP shunt
General surgery	6661EZE	Framed delivery incise drape		26cm x 20cm	Small hernia, appendectomy, port-a-cath placement (can be applied by single user)
	6648EU	Large incise drape		56cm x 60cm	Abdominal procedures
	6648EZE	Large incise drape		60cm x 60cm	Abdominal procedures
	6651EU	Extra-large incise drape		56cm x 85cm	Abdominal procedures
	6651EZE	Extra-large incise drape		60cm x 85cm	Abdominal procedures
	Cardiovascular	6635	Incise drape		10cm x 20cm
6661EZE		Framed delivery incise drape		26cm x 20cm	Fistula placement, clot removal, carotid (can be applied by single user)
6648EU		Large incise drape		56cm x 60cm	Lobectomy, valve replacement
6648EZE		Large incise drape		60cm x 60cm	Lobectomy, valve replacement
6651EU		Extra-large incise drape		56cm x 85cm	Abdominal procedures, AAA, fem-pop
6651EZE		Extra-large incise drape		60cm x 85cm	Abdominal procedures, AAA, fem-pop

*All dimensions are approximate.

Help reduce surgical site contamination risk.

Use this chart to see which 3M™ Ioban™ 2 Antimicrobial Incise Drape may be useful for your surgical application. Specialty drapes or general incise drapes may be appropriate based on product features and dimensions.

- Specialty antimicrobial incise drape for your application.
- Secondary option if a specialty drape is unavailable.

Procedure type		Drape type						Speciality incise drapes			
		Incise sheets									
		6635 Incise	6661EZE Framed delivery	6640EU/6640EZE Small	6650EU/6650EZE Medium	6648EU/6648EZE Large	6651EU/6651EZE XL	6617 Isolation	6619 Large isolation	6657 C-section pouch	6665 Abdominal-perineal
Orthopaedic	ACL										
	Arthroscopy										
	ORIF small extremity										
	Fractured neck of femur, hemi-arthoplasty, total hip replacement, intramedullary nailing										
	Paediatrics										
	Podiatry (including ankle)										
	Spine										
	Total hip										
	Total hip-anterior										
	Total knee										
	Total shoulder										
OB GYN	Open abdominal/pelvic surgery										
	C-section										
General	Abdominal-perineal resection										
	Appendicectomy										
	Colon resection (open)										
	Hernia repair										
	Open abdominal procedure										
	Liver transplant										
	Kidney transplant										
Cardiovascular/Thoracic	AAA (abdominal aortic aneurysm)										
	AV fistula										
	Coronary artery bypass										
	Carotid endarterectomy										
	Embolectomy										
	Femoropopliteal bypass										
	Lobectomy										
	Valve replacement										
	Pacemaker										
Neuro	Craniotomy										
	VP shunt										
	Spinal surgery										

For more information, contact your Solventum representative.



KCI Medical Limited
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB

3M Healthcare EMEA Export Limited
2 Cumberland Place
Fenian Street
Dublin 2
D02 H0V5

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.

© Solventum 2024. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates. 3M, the 3M logo and other 3M owned and/or licensed trademarks are trademarks of 3M Company. Other trademarks are the property of their respective owners. OMG1382784.